



METAL SUCTION TIPS

- Single use, for surgery and examination of the ear
- Sterile, single packed in peel pack
- ETO sterilised
- Special angle on the tip for ease of sight during examination and operations with microscope
- Quick connector system – it is easy to change tip during ongoing examination or surgery



<u>Art No</u>	<u>Size mm</u>	<u>Ward box / Trsp box</u>	
6066500042	0,4 x 60	lilac	60/900
6066500075	0,7 x 60	black	60/900
6066500109	1,0 x 60	brown	60/900
6066500141	1,4 x 60	white	60/900
6066500040	0,4 x 70	lilac	60/900
6066500070	0,7 x 70	black	60/900
6066500100	1,0 x 70	brown	60/900
6066500140	1,4 x 70	white	60/900
6066500208	2,0 x 70	green	60/900
6066500045	0,4 x 80	lilac	60/900
6066500072	0,7 x 80	black	60/900
6066500200	2,0 x 80	green	60/900
6066500250	2,5 x 80	orange	60/900
6066500300	3,0 x 80	blue	30/450
6066500400	4,0 x 80	yellow	30/450
6066509100	1,0 x 100	brown	30/900
6066509140	1,4 x 100	white	30/900
6066509200	2,0 x 100	green	30/900
6066509250	2,5 x 100	orange	30/900
6066509300	3,0 x 100	blue	30/900

ward box = min qty / order

Material: PP/ Medical grade
stainless steel ("virgin" steel)
MDD Rule: Rule 6
MDD Class: IIa
CE: Acc. existing directive

