

MEDI PLAST®



MEDIPLAST® NASAL SUCTION TIP WITH TRANSPARENT HANDLE

- For surgical operations and rhinological examinations
- Single use, sterile, single packed in peel-pack
- Medical grade stainless steel, transparent handle (PMMA)

Art.no.	mm	ward.box/trsp.box
66516120	2,0 x 160 mm	25/375
66516130	3,0 x 160 mm	25/375
66516140	4,0 x 160 mm	25/375

ward.box = min.qty/order

Material: PMMA/ Medical Grade
Stainless Steel (virgin steel)
MDD Rule: Rule 6
MDD Class: IIa
CE: Acc to regulation